



IRiS
Soft Switch
Instruction Manual

Congratulations on the purchase of your IRiS Soft Switch!

The IRiS Soft Switch has been designed specifically for people with special needs and should, if used in accordance with these guidelines, offer years of switching learning and fun!

The IRiS Soft Switch offers endless opportunities to develop numerous skills including colour recognition, cause and effect, hand to eye co-ordination, gross and fine motor skills and many more.

The IRiS Soft Switch is an IRiS Talker and can control any number of IRiS Listeners.

Put simply, it is operated by pressing any of the switches.

Please ensure you read this Instruction Manual in full before operating the unit.

Features

- IRiS Wirefree Compatible
- 8 Brightly Coloured Switches
- 6.5mm Switch Input Sockets
- 6.5mm to 3.5mm Jack Adaptors
- Rugged Switches
- 4 x AA Batteries
- Easy to use
- 2 Year Manufacturers Guarantee (Collect and Return)

Getting Started

1. Simply press any switch on the IRiS Soft Switch – the unit will control the products it is paired to.
2. To pair the IRiS Soft Switch to a product please see the 'Pairing to IRiS Listeners' section in this manual.

Using your own switches

Connect your own switches by plugging them into the sockets at the rear of the unit. The sockets are located in the same order as the colours on the top of the unit. 8 6.5mm to 3.5mm jack adaptors are supplied with the unit to ensure you can connect your own switch.

Pairing to IRiS Listeners

This procedure is split into 2 sections, one for those who do not have an IRiS Master Pairer and for those who do.

Without a IRiS Master Pairer

Every IRiS Listener has a button on it called a pairing button. To get the IRiS Soft Switch to work with a product, first decide which product you want to control then find the pairing button on that product and press it. You then need to press any switch on IRiS Soft Switch within 5 seconds. The IRiS Listener will now change colour every time you press a switch if correctly paired. The unit will now work with the IRiS Soft Switch.

With a IRiS Master Pairer

To get the IRiS Soft Switch to work with a product, first decide which unit you want to control and find its pairing button on the IRiS Master Pairer and press it. You then need to press any switch on IRiS Soft Switch within 5 seconds. The IRiS Listener will now change colour every time you press a switch if correctly paired. The unit will now work with the IRiS Soft Switch.

Cleaning and Maintenance

The unit should be cleaned with a damp soft cloth. A non abrasive cleaning agent may be used for stubborn stains. Ensure no liquid enters the unit.

The unit has been designed for use by people with special needs however this does NOT mean it is unbreakable. You must ensure the unit is not dropped, thrown, stood upon, immersed in water, etc.

Troubleshooting

Q. The unit doesn't work.

A1. When a switch is pressed on the unit does a blue LED illuminate? If not,

A2. Replace the batteries

A3. If the unit does not operate correctly after completing the above please contact us.

Q. The unit won't control any product.

A1. When a switch is pressed on the unit does a blue LED illuminate? If not,

A2. Replace the batteries

A3. Is the unit paired to an IRiS Listener? If so,

A4. Try pairing the unit to another IRiS Listener, does this work? If not, please contact us.

For support please email support@experia-innovations.co.uk or telephone 0845 644 0977. Please note: email is our support teams preferred communication and this takes precedence over telephone calls.